



# MEDPOR® Biomaterial Chin & Mandible Implants


Porex Surgical's Chin & Mandible Product lines made with MEDPOR Biomaterial

## CONTOURED TWO-PIECE CHIN IMPLANTS




CAT#	DESCRIPTION	A	B	C
86000	Contoured Two-Piece Chin	72mm x 42mm	x 3mm	
86001	Contoured Two-Piece Chin	74mm x 42mm	x 5mm	
86002	Contoured Two-Piece Chin	78mm x 50mm	x 7mm	
86003	Contoured Two-Piece Chin	80mm x 55mm	x 9mm	
85001	Chin Sizer Set for Contoured Two-Piece (Silicone, Non-Sterile)			

## TWO-PIECE CHIN IMPLANTS





CAT#	DESCRIPTION	A	B	C
8320	Small Projection	56mm x 33mm	x 5mm	
8321	Medium Projection	56mm x 36mm	x 7mm	
8322	Large Projection	57mm x 38mm	x 9mm	
9953	Chin Sizer Set for Two-Piece Design (Silicone, Non-Sterile)			

## RZ EXTENDED CHIN IMPLANTS



3mm Small  
5mm Medium  
7mm Large




CAT#	DESCRIPTION	A	B	C
8313	RZ Ext Round Chin - Small	45mm x 47mm	x 3mm	
8314	RZ Ext Round Chin - Med	45mm x 47mm	x 5mm	
8315	RZ Ext Round Chin - Large	45mm x 47mm	x 7mm	
8316	RZ Ext Square Chin - Small	45mm x 47mm	x 3mm	
8317	RZ Ext Square Chin - Med	45mm x 47mm	x 5mm	
8318	RZ Ext Square Chin - Large	45mm x 47mm	x 7mm	
9954	Chin Sizer Set for Extended Designs (Silicone, Non-Sterile)			


## BUTTON CHIN



Small (inset image)




Medium - Large




CAT#	DESCRIPTION	A	B	C
86010	Button Chin - Small	40mm x 25mm	x 4mm	
86011	Button Chin - Medium	47.5mm x 21.5mm	x 5.5mm	
86012	Button Chin - Large	48.5mm x 21.5mm	x 7mm	




## RZ MANDIBULAR ANGLE IMPLANTS



CAT#	DESCRIPTION	A	B	C
9955	Mandibular Angle RZ Left - Small	65mm	35mm	3mm
9956	Mandibular Angle RZ Right - Small	65mm	35mm	3mm
9957	Mandibular Angle RZ Left - Med	65mm	35mm	7mm
9958	Mandibular Angle RZ Right - Med	65mm	35mm	7mm
9959	Mandibular Angle RZ Left - Large	65mm	35mm	11mm
9960	Mandibular Angle RZ Right - Large	65mm	35mm	11mm
9566	RZ Mandibular Sizer Set (Silicone, Non-Sterile)			




## ANGLE OF THE MANDIBLE IMPLANTS




CAT#	DESCRIPTION	A	B	C	D	IP	LP
7537	Ramus w/Infer. Ridge E-5 - Left	79	32	5	10	5	7
7538	Ramus w/Infer. Ridge E-5 - Right	79	32	5	10	5	7
7539	Ramus w/Infer. Ridge E-10 - Left	79	32	10	10	10	7
7540	Ramus w/Infer. Ridge E-10 - Right	79	32	10	10	10	7
7541	Ramus w/Infer. Ridge ES-5 - Left	79	32	5	4	5	5
7542	Ramus w/Infer. Ridge ES-5 - Right	79	32	5	4	5	5
7543	Ramus w/Infer. Ridge ES-10 - Left	79	32	10	4	10	5
7544	Ramus w/Infer. Ridge ES-10 - Right	79	32	10	4	10	5
88051	Mand. Angle Ramus w/Infer. Ridge Lrg - Left	79	37	5	10	5	4
88052	Mand. Angle Ramus w/Infer. Ridge Lrg - Rgt	79	37	5	10	5	4

## CONTOURED MANDIBULAR ANGLE IMPLANTS



CAT#	DESCRIPTION	A	B	C	D
88037	Cont'd Mand. Angle - Left	59mm	29mm	7mm	11mm
88038	Cont'd Mand. Angle - Right	59mm	29mm	7mm	11mm

## LATERAL AUGMENTATION ONLY SHAPE



CAT#	DESCRIPTION	A	B	C	D	E	F
7535	Mandible Angle - Left	47	37	3	6.5	3	3
7536	Mandible Angle - Right	47	37	3	6.5	3	3
88053	Mandible Angle Large - Left	57	40	4	10	4	3
88054	Mandible Angle Large - Right	57	40	4	10	4	3

**POREX**  
SURGICAL PRODUCTS GROUP

MEDPOR is a registered trademark of Porex Surgical Inc.  
The Porex Surgical service mark is owned by or used under the license and authority of Porex Corporation.

SPG-1742-082807-01

Porex Surgical, Inc.  
15 Dart Road, Newnan, GA 30265-1017 USA  
1-678-479-1610 • FAX 1-678-423-1435

Porex Surgical GmbH  
Lindwurmstr 101, 80337 München Germany  
+49 (0)89 232415-0 • FAX +49 (0)89-23241515  
www.porexurgical.de

1-800-521-7321 • www.porexurgical.com

Please contact Porex Surgical, Inc.  
for a complete list of aesthetic and  
reconstructive MEDPOR Shapes.  
E-mail: surgical.info@porex.com