

POREX

SURGICAL PRODUCTS GROUP


Natural Iris Conformers

For Immediate Post-Op Enucleation Patients



Natural Iris Conformers

Natural Iris Conformers are a new rehabilitative conformer that allows the immediate post-operative enucleation patient to heal comfortably without an eye patch covering the socket for the time period required prior to fitting with a custom-made prosthesis.



“The Natural Iris Conformer is an excellent product for immediate post-op enucleation patients. Patients find them to be very comfortable. The biggest difference I have seen is the patient’s ability to get on with social activities before getting their custom prosthesis. I have had many patients that would wear patches or dark glasses. Some would wear their hair over their face or refuse to leave the house until they got their custom prosthesis. With the more normal appearing Natural Iris Conformer, patients feel more at ease.”

John D. Ng, M.D.

Indications for Use

Natural Iris Conformers are temporary devices used after the enucleation of the eye to maintain the shape of the socket and prevent closure or adhesions during the healing process. This product is for temporary use for approximately six to eight weeks following enucleation and then should be replaced by a custom impression-fitted ocular prosthesis fabricated by a trained ocularist.

Natural Iris Conformers are available in six colors: blue, blue-gray, green hazel, brown hazel, medium brown, and dark brown. Each color is available in three sizes: small, medium and large. Approximately 98 percent of patients require a medium conformer.



Comparison

The Natural Iris Conformer is slightly larger than standard clear conformers, but not quite as thick as a custom prosthesis.



"Two of my own children didn't even notice until I said something."



"The Natural Iris Conformer presents an immediate emotional and psychological lift at the most critical of times for the patient and their family."

Doss K. Tannehill, BCO





POREX

SURGICAL PRODUCTS GROUP

www.porexurgical.com

Porex Surgical, Inc.

15 Dart Road, Newnan, GA 30265-1017 USA

Email: surgical.info@porex.com

Telephone: 1-678-479-1610 • 1-800-521-7321

Fax: 1-678-423-1435

Porex Surgical GmbH

Balanstrasse 55

81541 München, Germany

Email: surgical.info2@porexurgical.com

Telephone: +49 (89) 232415-12

Fax: +49 (89) 232415-15

Copyright 2010, Porex Surgical Inc. The Porex Surgical service mark is owned by or used under the license and authority of Porex Corporation.
SPG-2098-020510-01