



August 20, 2002

New MEDPOR® Chin Implant From Porex Surgical Inc.

Porex Surgical Inc. is pleased to announce the addition of a new contoured two piece chin implant to the line of MEDPOR Biomaterial Implants available for facial reconstruction and aesthetics. The new MEDPOR Contoured Two Piece Chin Implant, designed with Michael J. Yaremchuck, M.D., provides an excellent anatomical fit to the bony anatomy with a gradual taper and concave posterior surface. The new chin implant is available in four sizes and enhances anterior projection at the mentum and subtle augmentation as it extends laterally along the ramus.

MEDPOR Implants are composed of biocompatible porous polyethylene. The interconnecting, omni-directional pore structure of MEDPOR Biomaterial allows integration with the patient's tissues by allowing fibrovascular ingrowth. MEDPOR Implants are available in over 250 standardized shapes, which can be modified by the surgeon for individual patient requirements.

For more information on the MEDPOR Contoured Two Piece Chin Implant or to request information about the complete line of MEDPOR facial implants, please contact Porex Surgical at 1-800-521-7321, or by e-mail at surgical.info@porex.com.

[Return to Press Release Index](#)