

February 18, 2002

## **New MEDPOR® Cranial Hemisphere Implant**

---

MEDPOR Cranial Hemispheres save time and cost in replacing bony loss of large cranial defects. This “off-the-shelf” alternative to customized implants, complex grafts, and other implant materials approximates the contour of the half cranium. The implants are 4mm thick and are available in left and right versions. They are provided sterile, individually packaged.

MEDPOR Implants are composed of biocompatible porous polyethylene. The interconnecting, omni directional pore structure of MEDPOR Biomaterial allows integration with the patient’s tissues by allowing fibrovascular ingrowth. MEDPOR Implants are available in over 250 standardized shapes, which can be modified by the surgeon for individually patient requirements.

For more information on the MEDPOR Cranial Hemisphere or to request a information about the complete line of cranial implants, please contact Porex Surgical at 1-800-521-7321, or by e-mail at [surgical.info@porex.com](mailto:surgical.info@porex.com).

[Return to Press Release Index](#)