

August 20, 2002

## **New MEDPOR® Orbital Rim Onlay Implant**

---

The new MEDPOR Orbital Rim Onlay Implant, from Porex Surgical Inc. was designed with Robert A. Goldberg, M.D. to offer support to the lower eyelid and lateral canthus in patients who have inadequate support due to congenital variation or changes associated with aging. The implant augments the inferior and lateral orbital rim and moderately increases the anterior rim projection. With the addition of the Orbital Rim Onlay Implant, Porex Surgical provides surgeons with a variety of styles to meet patient's individual requirements.

MEDPOR Implants are composed of biocompatible porous polyethylene. The interconnecting, omni-directional pore structure of MEDPOR Biomaterial allows integration with the patient's tissues by allowing fibrovascular ingrowth. MEDPOR Implants are available in over 250 standardized shapes, which can be modified by the surgeon for individual patient requirements.

For more information on the MEDPOR Orbital Rim Onlay Implant or to request information about the complete line of MEDPOR facial implants, please contact Porex Surgical at 1-800-521-7321, or by e-mail at [surgical.info@porex.com](mailto:surgical.info@porex.com).

[Return to Press Release Index](#)