



February 25, 2002

NovaBone-C/M™ Now Available in 1.5cc Packaging

Porex Surgical Inc. is pleased to announce the release of new smaller packaging for NovaBone-C/M, Synthetic Bone Graft Particulate. The new 1.5cc size is conveniently packaged in shelf packs of four individually sterile mixing cups to eliminate excess waste in the operating room when filling around bone grafts or burr hole areas.

NovaBone-C/M is a unique biocompatible synthetic bone graft particulate for use as a bone graft replacement or a bone graft extender in cases where particulate autograft or allograft might be used. The material is composed of Bioglass®, a radiolucent bioactive material consisting of sodium, calcium, phosphorous, and silicon oxides with the ability to bond to bone and soft tissue.

For more information on NovaBone-C/M or other Porex Surgical Inc. products please contact Customer Care at 1-800-521-7321, or by e-mail at surgical.info@porex.com.

[Return to Press Release Index](#)