

# MEDPOR® Biomaterial

For Oculoplastic Procedures



Additional resources available for your reference are:

- [www.porexurgical.com](http://www.porexurgical.com)
  - Product Information <http://www.porexurgical.com/english/surgical/skull.asp>
  - Clinical References <http://www.porexurgical.com/english/surgical/slitorb.asp>
- Surgical Videos <http://www.porexurgical.com/english/surgical/svideos.asp>
  - Enucleation with Use of the MEDPOR Conical Orbital Implant with Primary and Pre-placement of the Motility Coupling Post
  - Surgical Technique for Use of the MEDPOR Nasal Shell Implant
  - Natural Movement For Artificial Eyes
  - Orbital Floor Reconstruction Using MEDPOR® Biomaterial
- Suggestions for Use/Operative Techniques
  - BARRIER™ Sheet
  - SST™ Ocular Implants
  - MEDPOR® Coated Tear Drain
  - QUAD™ Motility Implant
  - COI Illustrative Guide
- Clinical Presentations CD in PowerPoint <http://www.porexurgical.com/english/surgical/clinicalpres.asp>
  - MEDPOR Properties and Handling Techniques
  - MEDPOR Motility coupling Post (MCP) Introduction
  - Lower Eyelid Spacer – Drs. C. Soparkar & J. Patrinely
  - Motility Coupling Post (MCP) – Secondary Placement – Drs. J. Choi and J. Shore

MEDPOR Surgical Implants are composed of porous high-density polyethylene. The porous structure of MEDPOR Implants supports ingrowth of soft and connective tissue. The mechanical properties of the MEDPOR material make it flexible and easy to reshape, yet sufficiently rigid for non-load bearing regions of the craniofacial skeleton. You may alter the shape of the material by cutting, carving, or relaxing the memory of the implant material by soaking it in heated sterile saline solution in the O.R.

Please contact Customer Care at 1-800-521-7321 to request additional information, non-sterile consultative samples or to place an order for your next procedure.

Thanks again for your interest.

---

**Porex Surgical, Inc.**

15 Dart Road, Newnan, GA 30265-1017 U.S.A. • 1-800-521-7321 • 1-678-479-1610 • FAX 1-678-423-1435

[www.porexurgical.com/surgical](http://www.porexurgical.com/surgical) • E-MAIL [surgicalinfo@porex.com](mailto:surgicalinfo@porex.com) • SPG-941-051002-01

The POREX SURGICAL service mark is owned by or used under the license and authority of Porex Corporation.