



**Newnan, Georgia – June 27, 2003**

## **New MEDPOR® Nasal Tip-Top™ Provides Predictable, Long Lasting Results**

---

The new MEDPOR Nasal Tip-Top Implant, from Porex Surgical Inc., is designed for patients who suffer from a flat or droopy nasal tip due to a congenital defect, trauma, or multiple rhinoplasties. The Nasal Tip-Top Implant allows the surgeon to provide predictable, long lasting results when native cartilage is missing or not available. The MEDPOR Nasal Tip-Top is uniquely designed to provide strength and flexibility while supporting and reshaping the nasal tip cartilages.

Designed with Dr. Paul O’Keeffe, the flat, wing shaped, 0.5mm thick implant features three strategically placed crimps for ease of shaping to create tip-defining points. After initial shaping the implant is placed over the tip cartilages to support or redefine the nasal tip. The degree of the angle of the folds determines whether the resulting tip is broad or more defined.

MEDPOR is a biocompatible porous polyethylene material. The interconnecting, omnidirectional pore structure allows for rapid fibrovascular ingrowth and integration of the patient’s tissue.

For more information on the MEDPOR Nasal Tip-Top Implant or to request an illustrative surgical guide please contact Porex Surgical Customer Care at 1-800-521-7321, or visit our website at [www.porexurgical.com](http://www.porexurgical.com).

###

**Contact:** Darlene Robinson  
Communications Manager  
Tel: 1-678-479-1610 Fax: 1-678-423-1435  
e-mail: [darlene.robinson@porex.com](mailto:darlene.robinson@porex.com)

#SPG-1193-062703-01