

MEDPOR® Biomaterial

For Reconstructive and Aesthetic Procedures



Additional resources available for your reference are:

- www.porexurgical.com
 - Product Information <http://www.porexurgical.com/english/surgical/skull.asp>
 - Clinical References <http://www.porexurgical.com/english/surgical/slitmidf.asp>
 - Surgical Videos on line <http://www.porexurgical.com/english/surgical/svideos.asp>
- Surgical Videos <http://www.porexurgical.com/english/surgical/svideos.asp>
 - MEDPOR RZ Mandibular Matrix System
 - Microtia Ear Reconstruction Using a MEDPOR Ear Implant
 - Microtia Reconstruction Using a MEDPOR Biomaterial
 - Nasal Dorsal Augmentation with MEDPOR Arch
 - Surgical Technique for Use of the MEDPOR Nasal Shell Implant
 - Cranioplasty Using MEDPOR Porous Polyethylene Implants
 - The MEDPOR Pterional Implant – A Neurological Approach
- Suggestions for Use/Operative Techniques
 - BARRIER™ Sheet
 - Reinisch Ear Procedure
 - External Nasal Valve Batten
 - Lateral Nasal Valve
 - RZ Mandibular Matrix System
- Clinical Presentations CD in PowerPoint <http://www.porexurgical.com/english/surgical/clinicalpres.asp>
 - MEDPOR Chin and Mandible Implants
 - MEDPOR Malar Implants
 - MEDPOR Custom Implants
 - MEDPOR Properties and Handling Techniques
 - Microtia Reconstruction – Dr. J. Reinisch
 - Temporalis Muscle Transfer Technique – Dr. Rapidis
 - Nasal Shell Implant – Dr. P. O’Keeffe
 - MEDPOR Properties and Handling Techniques

MEDPOR Surgical Implants are composed of porous high-density polyethylene. The porous structure of MEDPOR Implants supports ingrowth of soft and connective tissue. The mechanical properties of the MEDPOR material make it flexible and easy to reshape, yet sufficiently rigid for non-load bearing regions of the craniofacial skeleton. You may alter the shape of the material by cutting, carving, or relaxing the memory of the implant material by soaking it in heated sterile saline solution in the O.R.

Please contact Customer Care at 1-800-521-7321 to request additional information, non-sterile consultative samples or to place an order for your next procedure.

Thanks again for your interest.

Porex Surgical, Inc.

15 Dart Road, Newnan, GA 30265-1017 U.S.A. • 1-800-521-7321 • 1-678-479-1610 • FAX 1-678-423-1435

www.porexurgical.com/surgical • E-MAIL surgicalinfo@porex.com • SPG-943-051002-01

The POREX SURGICAL service mark is owned by or used under the license and authority of Porex Corporation.